



Patient Name / Essity File # _____ DOB _____

Address _____ Gender M F

City / State / Zip _____

Diagnosis _____ Date _____

Color	Quantity/Class	CCL 1 18-21 mmHg*	CCL 2 23-32 mmHg*	CCL 3 34-46 mmHg*
<input type="checkbox"/> Beige <input type="checkbox"/> Gray <input type="checkbox"/> Black <input type="checkbox"/> Cocoa <input type="checkbox"/> Cranberry <input type="checkbox"/> Cherry <input type="checkbox"/> Navy	Left			
	Right			
	Body Bandage CCL must be same as legs			

Styles

AD Knee AG-T Chap: pc. pr.
 AG Thigh AT Pantyhose
 B1G-T
 BG-T
 FT Biker Short

AT, BT, and B1T Pantyhose styles must be all one compression class. All leg lengths must be equal.

Straight Open Toe Length Slant Open Toe Length Slant Closed Toe Length
 Lateral _____ cm Medial _____ cm Medial _____ cm
 Straight Closed Toe Length Lateral _____ cm Lateral _____ cm
 Total Foot IZ _____ cm Total Foot IZ _____ cm

Options

T-Heel Top Comfort Zone (available for AG Thigh only)

Circum. (c)	Length (l)	Length (l)
cT ⁰	/GT	/T
cH ⁺⁺	/KT	/H

Tensions ** If measuring is done in lying position, cA please apply 0 tension


0 no tension
 + light tension
 ++ heavy tension

***cD/cG 0 tension with silicone band and straight ending

AT Panty Options

Adjustable Waistband
 Open Pubis
 Ribbed Fleece Waistband

Circumference (c)		Length (l): Taken from each landmark to floor	
Left	Right	Left	Right
		/K	
cG ⁺⁺		/G	
cF ⁺⁺		/F	
cE ⁺		/E	
cD ⁺⁺⁺		/D	
cC ⁺⁺		/C	
cB1 ⁺⁺		/B1	
cB ⁺		/B	
cY ⁰		/A (medial)	
cA ⁺		/A (lateral)	



Lateral Rise

AD standard 4cm
 Other: ____ cm (2cm-6cm)
 (AG fixed 7cm, no modifications)

Silicone Dot Top Band

2.5cm (AD Only) 5cm
 (Silicone band not available for AG-T)

Micro-Dot Top Band

5cm

SoftFit Band

5cm (AD Only)

Pocket	Lining (Pocket all sides closed)
<input type="checkbox"/> In-step	<input type="checkbox"/> In-step
<input type="checkbox"/> Back of knee	<input type="checkbox"/> Back of knee
	<input type="checkbox"/> Heel

All measurements should be in centimeters.

* Design Pressure

Comments:

Arion Easy-Slide

- The user-friendly donning aid for open toe compression stockings and tights
- The lightweight, smooth material provides ease of donning/application

Size	Shoe size	BNR	UOM / Box	Order Qty
X-Small	≤ 2	7965803	1	
Small	2.5 - 5.5	7965804	1	
Medium	6 - 8	7965802	1	
Large	8.5 - 11	7965902	1	
X-Large	≥ 11.5	7966001	1	



Arion Magnide®

- The ultimate donning aid for closed toe compression stockings and tights
- Unique magnetic closure for effective donning
- Has flexible sizing to fit different limb sizes

Size	Shoe size	BNR	UOM / Box	Order Qty
Medium	≤ 7.5	7957100	1	
Large	8 - 12	7957101	1	
X-Large*	≥ 13	7957102	1	



Arion Magnide® 2in1

- New and improved application aid
- For putting on and taking off closed toe compression stockings and tights
- All-in-one product
- User-friendly donning mechanism using own body weight

Size	Shoe size	BNR	UOM / Box	Order Qty
Medium	< 7.5	7957110	1	
Large	≥ 8 - 10	7957111	1	
X-Large*	≥ 10.5	7957112	1	



* XL also suitable for extra wide feet with cA up to 28 cm

Arion Sim-Slide®

- Multi-functional application aid
- For putting on and taking off open toe compression stockings and tights
- All-in-one product

Size	Shoe size	BNR	UOM / Box	Order Qty
Small	< 3 - 6	7063301	1	
Medium	6.5 - 8	7063302	1	
Large	8.5 - 10.5	7063303	1	
X-Large	≥ 11	7063304	1	



Arion Easy-Slide Kids

- Specifically designed to alleviate the problems experienced when donning arm and leg sleeves for children

Size	Age	Length	BNR	UOM / Box	Order Qty
X-Small	0-3 yrs	11.6" (29.5cm)	7957108	1	
Small	3-6 yrs	15.5" (39.5cm)	7957109	1	




Arion SlideX®

SlideX is especially helpful for patients with:

- Highly sensitive skin on the lower legs
- Open wounds with a plaster or bandage
- An above-average sized lower leg, caused by lipedema or lymphedema
- To be used in combination with Easy-Slide, Sim-Slide® and Magnide®

Size	Circumference of lower leg	BNR	UOM / Box	Order Qty
Medium	≤ 15.7" (40cm)	7957106	1	
Large	≤ 23.6" (60cm)	7957107	1	
Extra large	≤ 39.5" (100cm)	7957107	1	



Dycem® Non-Slip Material*

- Improves grip
- Provides stability
- Helps with easier donning of a JOBST® garment

* Caution: This product contains natural rubber latex which may cause allergic reactions.

Size	BNR	UOM / Box	Order Qty
One size	7763600	1	



Dycem® is a registered trademark of Dycem Ltd.

JOBST® Donning Gloves*

- Help to create friction and provide grip
- Help to evenly distribute the fabric on the leg

* Caution: This product contains natural rubber latex which may cause allergic reactions.

Size	BNR	UOM / Box	Order Qty
X-Small	7740027	1 pair	
Small	7740028	1 pair	
Medium	7740029	1 pair	
Large	7740030	1 pair	
X-Large	7740031	1 pair	



JOBST® Stocking Donner

- Easy-to-use device to help apply JOBST® Hosiery
- Ideal for people who have limited dexterity or have difficulty reaching their feet
- Use with closed and open toe styles
- Available in 2 sizes: Standard and Large calf

Size	BNR	UOM / Box	Order Qty
Standard	110913	1 donner	
Large	110910	1 donner	

