
sigvaris

Patient Name:		Contact Name:
Account Name:		Account #:
Account Phone #:		P.O. #:
Item #: 2739-AR	Height:	Weight:
Ship Name & Address:		Date:
Product Information		
Product includes one Chipsleeve Custom Arm, one pair of Cotton Liners, and one black Oversleeve.		
Right Arm		Color: Black



Important

Color: Black

Exact measurements are critical for this garment to ensure a proper fit. If you'd like to learn more about measuring and fitting, attend our MCE Basic Fitter Training. To sign up, visit **sigvariseducation.com** or contact your local territory manager for more information.

Alternatively, call Customer Care Solution Center at **800-322-7744**, or e-mail **us_orders@sigvaris.com**, to receive a remote consultation/training.

Supplies Needed

· Cell phone with camera. photos of the arm with measurement markings must be emailed to:

☐ Left Arm

- us_orders@sigvaris.com
- Measuring instructions and forms.
- SIGVARIS GROUP Measuring tape and body pen (or eyeliner pencil).
- · Signed Custom Order Terms & Conditions Form (include with order).

MEASURING INSTRUCTIONS

With patient seated, place the arm extended and elbow bent slightly, with the palm down, on a flat surface.

CIRCUMFERENCES

Measure circumference at wrist, and record on line A.

Mark dorsal aspect at distal edge of tape.

Measure circumference at forearm, and record on line A1.

Measure circumference at elbow, and record on line B.

Measure circumference at palm and record on line E.

Measure circumference at mid bicep, and record on line B1.

Measure circumference at axilla, and record on line C.

Mark dorsal aspect at proximal edge of tape.

LENGTH

Measure length at dorsal aspect from mark at Point A to mark at Point C and record in box D.







