IRC reduction kit^m order form - lower extremity


Customer Name



Insurance Covered Phone 309-664-7930 COMPRESSION Fax 309-664-7931 Email orders@medsourceorders.com



Customer No.	—— Purchase Order No. ———	Order Date	
Patient Name	Fax	Telephone	
Billing Address			
Shipping Address		Shipping method	
	- 11		
Measured by	Email	CC Info	



Le	g Measurement Guide		
G	Groin	В	Narrowest part of ankle above
E2	15cm above E		malleolus
E1	5cm above E	A1	Middle of foot
E	Center of patella	A	Ball of foot
D1	Knee crease	mall	Center of malleolus
D	Slightly below knee	1	ℓ = length
С	Widest part of calf	1	c = circumference

reduction kit lobe*: Indicate quantity in appropriate check boxes. * lobe is sold separately unless purchasing the reduction kit whole leg with lobe kit.

universal size additional ancillary items: x-wide undersleeve reduction kit leg (Pairs) shelf strap (Each) reduction kit toe cap*: Indicate quantity in appropriate check boxes. * toe caps are not included with kits and must be ordered separately.

	Short	Medium	Large	X-L	XX-L
сA	20-22 cm	22-24 cm	24-26 cm	26-28 cm	28-30 cm
Left					
Right					
Color	Beige		Black		

whole leg kit: Check knee or lobe kit option, and indicate corresponding length & width.

knee

reduction kit lobe

lower leg: Check the box next to the corresponding length and width. *cC cannot be more than 30cm larger than cB.

		Wi	dth				Lengt	h			
	Regular		Wide		Short		Standard		Long		Quantity
Left	℃ < 65 cm		℃ < 80 cm		ℓ mall-D1		ℓ mall-D1		ℓ mall-D1		
Right	CB < 55 cm		cB< 65 cm		< 32 cm		32 - 36 cm		> 36 cm		

upper leg: Check the box next to the corresponding length and width. *cG cannot be more than 30cm larger than cE2.

		Wi	dth				Length						
	Regular		Wide		Short		Standard		Long		Quantity		
Left Right	¢G < 80 cm ¢E2 < 65 cm		℃G < 125 cm ℃E2 < 110 cm		ℓ D1-G < 37 cm		ℓD1-G 37-41 cm		ℓ D1-G > 41 cm				
Each up leg kit in			 1 lower leg 1 knee or lo 1 upper leg 	, bbe co	mponent	who	ir of under ble leg extr ielf straps			hook s BTS ca			

u *	pper cG cann	leg kit	t: Che e than	eck the box 30cm larg	next ger tha	to the corr an cE2.	espor	nding length	and w	idth.		Each upper leg kit includes: • 1 reduction upper leg
			component									
	Short Sta		Standa	Standard Long			Regular Wide				 1 reduction shelf strap 6 hook and loop stays 	
	Left	ℓ D1-G		ℓ D1-G		ℓ D1-G		¢G < 80 cm		¢G < 125 cm		 1 Built-In-Tension system guide card
	Right	< 37 cm		37- 41cm		> 41 cm		¢E2 < 65 cm		<i>C</i> E2 < 110 cm		 1 paper measuring tape 1 direction for use

knee kit: (the reduction knee kit is only available in one length & width)*cE2 cannot be more than 30cm larger than cC.

	Length	Width	Quantity	Each knee kit includes:
Left	30 cm	¢E 2 < 110 cm ¢E < 102 cm		 1 reduction knee comp 2 hook and loop knee s
Right	50 GM	CC < 80 cm		• 1 paper measuring tap

e component • 1 reduction shelf strap knee spines ring tape

• 1 direction for use

lower leg kit: Check the box next to the corresponding length and width. *cC cannot be more than 30cm larger than cB.

			Leng	th				Wi	dth		Each lower leg kit includes:
	Shor	t	Standa	ard	Long	3	Regula	r	Wide	Quantity	• 1 reduction lov
Left	ℓ mall-D1		ℓ mall-D1		ℓ mall-D1		℃ < 65 cm		¢C< 80 cm		component • 1 pair reductio
Right	< 32 cm		32 - 36 cm		> 36 cm		¢B< 55 cm		cB< 65 cm		undersleeve leg • 1 reduction she

foot options: Indicate quantity in appropriate check boxes.

foot opt	ions are not ir	nclude	ed with kit	san	id mu	st be o	rdered	separa	itel	у.				
	customiza interlocki afw			npre ankl	essio let	n		ompre inders			pac band™	,	single band afv	N
	One Siz	e	Standa	rd	La	rge	Standard			Large	X-Large	e	One Size	e
¢ A1	< 44 cm		19-36 ci	m	22-4	l1 cm	19-3	2 cm	24	4-41 cm	< 44 cm			
с В	< 56 cm	< 56 cm		m	25-42 cm		19-2	7 cm	25	5-42 cm	< 56 cm			
Left														
Right														
reduction kit ext		exte	ension	:							foam pa	ad	lding:	
Shor (30cr			ndard 5cm)			Lo (40	ong Icm)				flec pad (2 pk)			t

werleg on kit


- g
- elf strap
- 6 hook and loop stays • 1 Built-In-Tension
- system guide card
- 1 paper measuring tape


• 1 direction for use

- tec pad (2 pk)

IRCAID reduction kit™ order form - upper extremity

Customer Name			mear
Customer No	Purchase Order No	Order Date	
Patient Name	Fax	Telephone	
Billing Address			
Shipping Address		Shipping method	
Measured by	Email	CC Info	





A	m Measurement Guide
G	Axilla
F	Midpoint of E and G
Е	Elbow crease
D	Midpoint of C and E
С	Wrist
В	Base of thumb webbing
	ℓ = length
	c = circumference

arm kit: Check the box next to the corresponding length and width. *If cB is > 30cm hand wrap will not t, though reduction arm component may still be used.

		Ler	ngth								
	Standa	ard	Long	Long		Regular		Wide		X-Wide	
Left Right	ℓ C-G < 45 cm		ℓ C-G > 45 cm		€ F< 50 cm € B< 30 cm*		CF 50-65 cm CB< 30 cm*		¢F 65-80 cm ¢B< 30 cm*		

Each upper arm kit includes:

1 reduction arm component
1 pair undersleeve arm

· 1 customizable hand wrap

1 Built-In-Tension guide card
1 direction for use

additional ancillary items:

undersleeve arm standard	undersleeve arm wide	reduction kit	
Standard <50cm max circ.	Wide <80cm max circ.	shelf strap	
Quantity (Each)	Quantity (Each)	Quantity (Each)	

reduction kit glove*:

Indicate quantity in appropriate check boxes. *not included in arm kit

	Short	Medium	Large	X-L	XX-L
С В	18-22 cm	20-23 cm	22-25 cm	24-27 cm	26-29 cm
Left					
Right					
Color	Beige		Black		

reduction kit vest: Check the box next to the corresponding length and width. Note, width has two measurement criteria. Select the size that meets both criteria.

	Width				Length (Axilla-Waist)						
	Sma	II	Larg	e	Larg	e	Standa	ard	Long	ł	Quantity
Chest Circum.	70 - 120 cm		121 - 180 cm		80 - 180 cm		20 - 30 cm		34 - 47 cm		
Chest Width (Axilla-Axilla Contour)	≤ 39 cm		35 - 39 cm		> 39 cm		20 - 30 cm				

Each vest kit includes:

• 1 reduction vest component 1 flec (at localized edema control) foam pad

1 direction for use

foam padding:

	flec pad (2 pk)		tec pad (2 pk)	
--	--------------------	--	-------------------	--

reduction head & neck:

universal size

Each head & neck kit includes:

• 1 large neck strap 1 tec (textured edema control) foam pad

• 1 small facial chin strap 1 direction for use 1 direction for use

P1515 Rev